

# Summer Internship : Brain PET Image Denoising

2022.08.19

Soyon Park



SEOUL  
NATIONAL  
UNIVERSITY

Functional & Molecular  
Imaging System Lab



# Outline

## Brain PET Image Denoising

1. Introduction
2. Method
3. Results
4. Discussion
5. Conclusion
6. Future work

# 1. Introduction


- **PET (Positron Emission Tomography) Imaging**
  - **Low spatial resolution** and **signal-to-noise ratio** limit the **detection** and **quantitative accuracy** of PET imaging



⇒ **To improve PET image quality, denoising** has become a **crucial pre-processing step** for PET imaging

# 1. Introduction

- **PET image:** three-dimensional, volumetric medical data
- **3D Medical Image**




# 1. Introduction

- Learning Algorithm

## Noise2Clean

$$\operatorname{argmin}_{\theta} \sum_i L(f_{\theta}(\hat{x}_i), y_i)$$

$f_{\theta}$  = network function,  $\hat{x}_i$  = noisy input,  $y_i$  = clean target




Supervised Learning

## Noise2Noise

$$\operatorname{argmin}_{\theta} \sum_i L(f_{\theta}(\hat{x}_i), \hat{y}_i)$$

$f_{\theta}$  = network function,  $\hat{x}_i$  = noisy input,  $\hat{y}_i$  = noisy target




Self-supervised Learning

# 2. Method

## 1. Dataset

- **Dataset:**  $^{18}\text{F}$ -Florbetaben Brain PET data of 47 patients (40 for Training/ 7 for Testing)
- $\beta$ -amyloid neuritic plaque density in the brain  $\Rightarrow$  Alzheimer's disease diagnosis
- $^{18}\text{F}$ -Florbetaben is a radiotracer highly specific and sensitive for the  $\beta$ -amyloid neuritic plaque density
- **The intake of  $\beta$ -amyloid deposits in  $^{18}\text{F}$ -FBB informs that it is the progress of Alzheimer's disease**




Normal and positive  $^{18}\text{F}$ -florbetaben scans in regions of interest.

# 2. Method


## 1. Dataset

- List-mode PET data  
: 20 data bins with 30 second duration and 600 second reference data
- PET image matrix size:  $200 \times 200 \times 109$  , pixel size:  $2.04 \times 2.04 \times 2.03 \text{ mm}^3$



# 2. Method

## 2. Network Architecture and Training







# 2. Method

## 2. Network Architecture and Training

### 2D U-Net




- 2D Data Preprocessing



# 2. Method

## 2. Network Architecture and Training

### 2D U-Net




# 2. Method

## 2. Network Architecture and Training

### 2D U-Net

④ Transaxial+Coronal+Sagittal




# 2. Method

## 2. Network Architecture and Training

### 2D U-Net


⑤ Tri-planar (3 channel)



# 2. Method

## 2. Network Architecture and Training

### 3D U-Net



- **3D Data Preprocessing**

- $32 \times 32 \times 32$  patch
- A total of 3200 possible patch center positions are calculated by applying a small random perturbation to the equally distanced 3D grids
  - Randomly select 400 patches
- Patch Thresholding:



$$p = (1 - |r - 0.5| \times r)$$

$p$ : The relative probability that a patch center position is selected,  $r$ : ratio of foreground pixels in the patch

# 2. Method

## 2. Network Architecture and Training

### ACS U-Net



# 2. Method

## 2. Network Architecture and Training

### Normalization

- 99<sup>th</sup> percentile

⇒ Normalize wide dynamic range of PET image intensity (0~thousands) → (0~1.xx)

### Training

- Batch size: 16
- Epoch: 30
- Learning rate: 0.0001(Initial)/ “ReduceOnPlateau” learning rate scheduler
- Optimization: Adam
- Clinical PET images: L2 loss (mean square error)

## 2. Method

- **Models**
- **2D/3D/ACS Convolution**
  - **2D U-Net**
    - Transaxial
    - Coronal
    - Sagittal
    - Transaxial+Coronal+Sagittal
    - Tri-planar: (Transaxial, Coronal, Sagittal) 3 channels
  - **3D U-Net**
    - 3D patch based
  - **ACS U-Net**
    - 3D patch based
- **Learning Algorithm**
  - **Noise2Clean**
  - **Noise2Noise**



## 2. Method

### 3. Image Analysis

#### 1) PSNR

$$PSNR = 10 \log_{10} \left( \frac{MAX^2}{MSE} \right)$$

*MAX: maximum value of the reference image, MSE: mean square error between the tested and reference image*

- **PSNR**: common method to **measure the fidelity** that is independent of the dynamic range of images

#### 2) SSIM

$$SSIM = \frac{(2\mu_x\mu_y + C_1)(2\sigma_{xy} + C_2)}{(\mu_x^2 + \mu_y^2 + C_1)(\sigma_x^2 + \sigma_y^2 + C_2)}$$

*x, y: compared images,  $\mu$ : mean,  $\sigma$ : variance,  $C_1, C_2$ : constants*

- **SSIM**: take advantage of the **human visual system**, assessing image quality by extracting **structural information** and **calculating the similarity** by three comparisons: **Luminance, Contrast, and Structure**

# 3. Results

## PSNR: 2D/3D U-Net N2C/N2N


U-Net	Data	Noise2Clean		Noise2Noise	
Input		33.11			
2D U-Net	Transaxial	36.85	37.20	36.57	37.20
	Coronal	37.74		37.86	
	Sagittal	37.01		37.17	
	Transaxial +Coronal +Sagittal	37.08		36.15	
	Tri-planar (3 channels)	36.17		36.28	
3D U-Net	3D Patch	37.36		37.34	
ACS U-Net	3D Patch	37.43		37.41	

# 3. Results


## SSIM: 2D/3D U-Net N2C/N2N

U-Net		Data	Noise2Clean		Noise2Noise	
Input			0.9893			
2D U-Net	Transaxial	0.9954	0.9956	0.9951	0.9955	
	Coronal	0.9959		0.9956		
	Sagittal	0.9954		0.9957		
	Transaxial +Coronal +Sagittal	0.9948		0.9948		
	Tri-planar (3 channels)	0.9946		0.9944		
3D U-Net	3D Patch	0.9972		0.9971		
ACS U-Net	3D Patch	0.9975		0.9974		

# N2C Results




# N2N Results



# 4. Discussion

- **Region of Interest (ROI)**

- Training the image only with the ROI → improve the performance of denoising
- Brain would be the ROI and background should be excluded
- If the background of an image is included when measuring PSNR and SSIM  
⇒ Quality of noisy image can be overestimated



# 5. Conclusion

- We compared **2D U-Net** and **3D U-Net** for **Noise2Clean** and **Noise2Noise** framework to reduce the noise in short scan time  $^{18}\text{F}$ -Florbetaben brain PET images
  - The applied deep learning methods showed **remarkable performance** for **reducing noise** in PET images
  - Noise2Noise gave comparable results to Noise2Clean without using reference images
- ⇒ **Self-supervised methods** can replace the traditional supervised method

# 6. Future Work

- **Other Convolution Methods**

- Transfer Learning
  - Hybrid Convolution: H-Dense U-Net
    - Uses 2D pretrained networks with multi-slice inputs and 3D networks that are initialized randomly with volumetric inputs together for training

- **Other Learning Methods**

- Self-supervised Methods
  - Noiser2Noise
    - Image denoising without access to clean training examples or access to paired noisy training examples
    - requires only a single noisy realization of each training example and a statistical model of the noise distribution



# Thank you for your attention

For more information,  
Codes on GitHub: <https://github.com/soyon007>  
Email Address: [soyon0304@naver.com](mailto:soyon0304@naver.com)